

# Fendall Eyesaline® Wall-Mount Eye Wash Station

PLS1497 Wall-Mount; Not Plumbed, Holds 32 oz., Flow Rate Variable

Economical stations can be mounted throughout your facility to soothe eyes with saline solution until you can reach a primary eye wash.

- Easy-to-clean station mounts directly to wall with included hardware; molded-in overhang protects bottle from dust and debris
- Saline solution is preserved and buffered to match natural pH chemistry of human tears
- Sterile, FDA-grade polyethylene bottle is sealed with a break-off T-tab lid for tamper resistance; easy to open in an emergency
- Extended-flow nozzle provides irrigation longer than conventional designs
- Elastic band secures bottle in place when not in use
- Hi-viz green station stands out in low lighting and when vision is blurred
- Long shelf life so you can keep a replacement in-stock following an emergency; clearly marked expiration date makes inspection fast and foolproof
- Features bilingual text on bottle and full English instructions on station



SPERIAN

## Specifications

Capacity	32 oz.
Style	Wall-Mount; Not Plumbed
Color	Green
Dimensions	3.25" W x 12.25" L x 14" H
Dispensing Method	Bottle
Fluid Storage Method	Bottle
Sold as	1 each
Weight	2.84 lbs.
Composition	Station - Plastic Solution - Saline
Flow Rate	Variable
Includes	1 - Wall-Mount Station

1 - Mounting Hardware

1 - 32 oz. Bottle of Solution (PLS1499)

**Shelf Life**

3 years

**UNSPSC**

46181810

**Pigalog® Page Number**

Page 224

## Metric Equivalent

**Dimensions**

8.3cm W x 31.1cm L x 35.6cm H

**Weight**

1.3 kg

## Technical Information

### CountryShipOnly

This product can ONLY be shipped to: United States, Australia, Iraq.

## Technical Documents

Sterile Eyesaline® Solution

29 CFR 1910.151(c)



**New Pig**

World's best stuff for leaks, drips and spills.

**One Pork Avenue • Tipton, PA 16684-0304**

1-855-493-4647 • Fax: 1-800-621-7447 • [newpig.com](http://newpig.com) • [hothogs@newpig.com](mailto:hothogs@newpig.com)